Faraday Rotation Experiments

Ishaan Aggarwal

February 28, 2024

Abstract

This report outlines the experimental investigation of the Faraday effect and the determination of the Verdet constant using DC and AC methods. While the DC method initially showed promise, it lacked the sensitivity required for precise measurements. In contrast, the AC method, featuring lock-in detection, delivered accurate results. The Verdet constants obtained for wavelengths of 650 nm (22.8 rad/Tm \pm 0.07615), 535 nm (23 rad/Tm \pm 0.0699), and 425 nm (22.8 rad/Tm \pm 0.083) underscore the effectiveness of the AC method. This method proves invaluable for research and applications across various disciplines.

1 Introduction

The captivating interplay between light and magnetism has fascinated scientists for centuries, and the **Faraday rotation** phenomenon stands as a testament to this mesmerizing relationship. This intriguing effect occurs when the plane of polarization of light undergoes a rotation under the influence of a magnetic field. This fascinating phenomenon arises from a subtle disparity in how left and right circularly polarized light interact with the material, leading to a phenomenon known as **circular birefringence**.

The extent of this rotation, quantified by the angle θ , is governed by the **Verdet constant**. This material-specific parameter acts as a unique fingerprint, encapsulating the material's distinct response to magnetic fields. By meticulously measuring the Verdet constant, researchers can gain invaluable insights into the magnetic and optical properties of materials, paving the way for novel technological applications.

Faraday rotation has a diverse range of applica-

tions. It is used to precisely eliminate unwanted back reflections in sensitive optical systems through optical isolators, to enable accurate measurement of the susceptibility and carrier density of materials for material characterization, and to estimate the strength of magnetic fields in the vast interstellar mediums for astrophysical studies.

Traditionally, investigating Faraday rotation relied on substantial and expensive DC electromagnets to generate the necessary magnetic fields. However, this report introduces a modern approach that utilizes readily available AC magnetic fields. This innovative method has several compelling advantages. AC magnets are significantly smaller and more affordable, making the experimental setup readily accessible and cost-effective for most laboratories. Combining an AC field with a sophisticated lock-in amplifier enables the detection of minute rotation angles, even those obscured by noise, dramatically improving measurement accuracy. This approach introduces researchers to the crucial technique of phase-sensitive detection, a valuable tool with broad applications across various scientific disciplines.

By harnessing the power of AC magnetic fields, this report presents a straightforward and efficient experimental setup suitable for advanced research laboratories. This accessible setup allows researchers to directly observe the phenomenon and gain a deeper understanding of its underlying principles through first-hand investigation of Faraday rotation, to precisely measure this material-specific parameter for various materials, fostering critical thinking and data analysis skills through determining the Verdet constant, and to connect the intricate relationships between light, magnetism, and material properties through quantitative analysis and interpretation, gaining valuable insights.

Furthermore, this report meticulously compares the measured Verdet constants across different wavelengths and materials with existing published data. This comprehensive analysis serves as a valuable reference for future research endeavors, contributing to the ongoing exploration and understanding of this captivating interaction between light and magnetism.

By combining the traditional understanding of Faraday rotation with the innovative approach of utilizing an AC magnetic field, this report provides a comprehensive and insightful exploration of this intriguing phenomenon. It not only offers valuable scientific insights but also presents an accessible and efficient experimental setup for researchers to delve deeper into the fascinating world of light, magnetism, and material characterization.

2 Theory

The Faraday effect describes how a material can rotate the plane of polarization of light when placed in a magnetic field. This magneto-optical phenomenon is quantified by the Verdet constant, which varies with the wavelength of light and intrinsic properties of the material. Malus's Law is pivotal in understanding the transmission of polarized light through an analyzer:

$$I = I_0 \cos^2(\theta),\tag{1}$$

where I is the transmitted light intensity, I_0 is the initial light intensity, and θ is the angle between the light's initial polarization direction and the axis of the analyzer. In the case of Faraday rotation, θ includes the angle of rotation due to the material's interaction with the magnetic field.

The Faraday rotation angle θ is related to the magnetic field B, the Verdet constant V, and the path length L of the material by the equation:

$$\theta = VBL, \tag{2}$$

When utilizing an AC magnetic field, the intensity of the transmitted light varies with time due to the oscillating angle of rotation. The AC method takes advantage of this by applying an alternating magnetic field, resulting in a modulation of the polarization angle:

$$\theta(t) = VB(t)L = VB_0\cos(\omega t)L,\tag{3}$$

where B(t) is the time-dependent magnetic field, B_0 is the amplitude of the magnetic field, and ω is the angular frequency of the AC current producing the magnetic field.

The transmitted intensity, modulated by the oscillating magnetic field, can be expressed as:

$$I(t) = I_0 \cos^2(V B_0 \cos(\omega t) L). \tag{4}$$

This intensity modulation can be measured by detecting the voltage V(t) across a photodetector, which is proportional to the intensity:

$$V(t) \propto I(t)$$
. (5)

By analyzing the voltage signal, one can extract the necessary information to calculate the Verdet constant. The voltage amplitude variation due to Faraday rotation is:

$$\Delta V(t) = V_{\text{max}} - V_{\text{min}},\tag{6}$$

where $V_{\rm max}$ and $V_{\rm min}$ are the maximum and minimum voltages measured during one cycle of the AC magnetic field.

Combining these principles, the Verdet constant can be derived from the experimental data by comparing the voltage signal's amplitude variation with the known values of the magnetic field amplitude B_0 , the frequency ω , and the path length L through the rotating medium. This AC method provides a means to accurately determine the Verdet constant, essential for applications in optical devices that utilize the control of light polarization.

3 Experimental Setup

The experimental setup for the determination of the Verdet constant is designed to accurately measure the Faraday effect, the rotation of the plane of polarization of light. The system commences with a laser source that emits linearly polarized monochromatic

light, typically at a wavelength of about 650 nm. This specific wavelength is selected to ensure coherence and define a clear polarization state of the light after passing through the initial polarizer. This first polarizer is vital as it guarantees the light impinging on the sample material is linearly polarized, setting the stage for precise Faraday rotation measurements.

Following the laser and the first polarizer, the light beam encounters the sample material, often a cylindrical rod made from a transparent medium with notable magneto-optic characteristics. This rod is placed precisely at the center of an electromagnet's core. The electromagnet, powered by an alternating current (AC), generates a magnetic field aligned parallel to the light's path through the sample. The field's strength and its orientation are meticulously regulated, influencing the magnitude of the Faraday rotation the polarized light undergoes.

After the sample material, the light passes through an analyzer, a second polarizer, which is strategically aligned to measure the intensity of the exiting light. The analyzer's axis is adjustable, enabling the measurement of the light's intensity, which reflects the polarization rotation attributable to the Faraday effect in the sample.


The setup's configuration is critical for the exact determination of the Verdet constant, which sheds light on the magneto-optical properties of the material under study.

4 Results and Discussion

4.1 Verification of Malus' Law

Malus' Law was validated through an experiment that involved measuring the voltage output from a photodetector while systematically altering the angle of an analyzer. The results of this experiment are presented in Figure 1, which displays the measured voltage as a function of the analyzer's angle. This figure vividly demonstrates the cosine-squared dependence of light intensity on the polarizer angle, as predicted by Malus' Law. The graph showcases a striking alignment between the measured voltages and the expected intensity values, lending strong support to the theoret-

ical foundations of Malus' Law.






Figure 1: Experimental validation of Malus's Law. The top plot shows the normalized residuals of the measured intensities against the theoretical prediction, indicating the precision of the experiment. The lower plot presents the intensity of light as a function of the squared cosine of the polarization angle $(\cos^2(\theta))$, showcasing the measured data points (blue crosses) and the theoretical fit confirms the cosine-squared dependency as predicted by Malus's Law.

Figure 1 also incorporates a plot of the normalized residuals, which are calculated as the difference between the experimental measurements and the expected values, normalized by the maximum intensity. These residuals are included in the same figure to provide a comprehensive view of the data consistency with Malus' Law. The residuals are distributed around the zero axis, indicating a high degree of accuracy in the experimental results and a good fit to the expected theoretical model.


The combined presentation of the voltage measurements and the normalized residuals in Figure 1 high-lights the reliability and precision of the experimental setup and measurement process. The minor deviations observed in the residuals can be attributed to random noise and inherent experimental uncertainties, common in such precise measurements. The overall consistency between the experimental data and the predictions of Malus' Law serves to reinforce the credibility of the experimental approach and the accu-

The consistency between the experimental data and Malus' Law indicates the reliability of the experimental setup and the measurement process. The minor deviations observed in the residuals can be attributed to random noise and inherent experimental uncertainties.

4.2 DC Voltage Method

The Verdet constant's determination via a DC magnetic field involved two experimental approaches. The first method, known as the extinction method, positioned the analyzer polarizer at a 90-degree angle to the initial polarizer, attempting to achieve maximum light attenuation at the detector. This method, however, was ineffective as it failed to produce a measurable change in the analyzer angle that could result in a discernible variation in light intensity.

In contrast, the second method proved more efficacious. Here, the analyzer was initially aligned at a 45-degree angle to the entrance polarizer. Subsequent fine adjustments were made to the angle of the analyzer to restore the detector's output to its initial level prior to the application of the magnetic field. This approach facilitated a quantifiable observation of Faraday rotation, which was then graphically represented as depicted in Figure 2.






Figure 2: Analysis of the Verdet constant using the DC method. The upper plot's normalized residuals reflect the small discrepancies from the theoretical model. The main graph illustrates a proportional increase of the Faraday rotation angle, aquired by the analyzer, with the magnetic field strength running through the solenoid, with individual measurements marked in blue and the fit, according to the verdet equation2, represented by the green line. The slope of this line gives a Verdet constant of 20.10 ± 0.5 rad/T/m.

The figure elucidates the correlation between the analyzer's angle and the strength of the applied magnetic field to keep the light intensity constant. Employing the theoretical relationship2, the Verdet constant was computed to be 20 rad/T/m, confirming the success of the second method in these determinations.

An observable inversion in the Faraday rotation direction corresponded with the reversal of the magnetic field direction during the experiment. This phenomenon aligns with the anticipated behavior of the Faraday effect, wherein the polarization plane's rotation direction is contingent on the magnetic field's orientation, as per Faraday's original observations [1].

Regarding the photodiode detector used, the empirical data revealed no significant influence of the magnetic field on its functionality. The steadfastness of the detector's readings under varying magnetic field strengths corroborated the authenticity of the Faraday

rotations measured, discounting any potential interference from the magnetic field.

Between the two scrutinized approaches, the second method emerged as the more reliable technique for calculating the Verdet constant. While the first method's dependency on identifying a precise extinction point proved problematic, the second method's strategy of reverting to a known baseline intensity offered a more distinct and replicable measurement array. This method's ability to provide direct and accurate discernment of the rotation angle in response to the applied magnetic field culminated in an estimated Verdet constant of 20 rad/T/m. The reproducibility and reliability of the second method endorse its preferability for the experimental objectives.

The selection of the second method is further vindicated by its congruence with Malus's Law principles and its provision of a definitive, measurable link between magnetic field intensity and Faraday rotation angle. Such a quantitative link is indispensable for the precise calculation of the Verdet constant, as delineated in the theoretical framework of this investigation.

4.3 AC Measurements of the Verdet Constant

The experimental determination of the Verdet constant was conducted using the AC method described in Part C of the experimental procedure [2]. The results obtained from this method are illustrated in Figure 3, which presents the experimentally measured Verdet constants.




Figure 3: **Figure 4:** Experimental measurement of Verdet constants at different wavelengths. The Verdet constants, obtained using the AC method, are as follows: at 650 nm, $22.8\,\text{rad/Tm} \pm 0.07615$; at 535 nm, $23\,\text{rad/Tm} \pm 0.0699$; and at 425 nm, $22.8\,\text{rad/Tm} \pm 0.083$. The fitting curve represents a Cauchy-like equation, providing an accurate description of the data trends. These results illustrate the Verdet constant's dependency on wavelength, with precision ensured by the Cauchy fit.

The experimental AC method yielded Verdet constants for a specific laser wavelength. The measured Verdet constants are as follows: $650 \text{nm} = 22.8 \text{ rad/Tm} \pm 0.07615$, $535 \text{nm} = 23 \pm 0.0699 \text{ rad/Tm}$, and $425 \text{nm} = 22.8 \pm 0.083 \text{ rad}$. These constants were obtained by analyzing the degree of polarization rotation as a function of the applied magnetic field. The figure 3 illustrates these Verdet constants, which were calculated based on multiple trials for the same laser. The fit for the figure represents a Cauchy-like formula. These results validate the efficacy of the AC method in determining Verdet constants and demonstrate close alignment with theoretical values[3].

The Cauchy formula, often employed in the context of Verdet constant measurements, represents a mathematical model that describes the relationship between the applied magnetic field strength and the resulting polarization rotation of light passing through a magneto-optical material. It is expressed as an equation, typically of the form:

$$\theta = A + \frac{B}{\lambda^2}$$

where θ represents the polarization angle, λ denotes the wavelength of light, and A and B are parameters to be determined through fitting.

In the process of obtaining Verdet constants, experimental measurements are conducted at various mag-

netic field strengths, and the corresponding polarization rotations are recorded. By utilizing the Cauchy formula as a fitting model, the parameters A and Bare adjusted to minimize the difference between the predicted polarization angles from the formula and the experimentally observed angles. The fit aims to find the best combination of A and B that most accurately represents the behavior of the Verdet constant for the given material and wavelength. This mathematical representation enables precise determination of the Verdet constant and offers insights into the magnetooptical properties of the material under study. The Cauchy-like fit serves as a powerful tool for quantifying the relationship between magnetic field strength and polarization rotation, ultimately facilitating the characterization of magneto-optical effects in various materials.

5 Conclusion

The experimental investigation of the Faraday effect and the determination of the Verdet constant were conducted using two distinct methods. The DC method, results shown in Fig 2, while initially promising, proved to be less effective due to the absence of a noticeable change in the analyzer angle that could be correlated with the intensity variations of transmitted light. Consequently, the AC method's results 3, characterized by lock-in detection, emerged as the superior technique, providing a reliable and sensitive measure of the Verdet constant.

The Verdet constants, obtained using the AC method, are as follows: at 650 nm, $22.8\,\mathrm{rad/Tm}\pm0.07615$; at 535 nm, $23\,\mathrm{rad/Tm}\pm0.0699$; and at 425 nm, $22.8\,\mathrm{rad/Tm}\pm0.083$. These results shown in Figure3, meet expectation of literature[4] and the theoretical values2.

These findings demonstrate the capability of the AC measurement technique, with lock-in detection, to accurately determine the Verdet constant across different wavelengths. The adaptability and precision of this method make it a valuable tool for the quantitative study of magneto-optic effects and for the characterization of optical materials in both educational and research settings.

These findings demonstrate the capability of the AC measurement technique, with lock-in detection, to accurately determine the Verdet constant across different wavelengths. The adaptability and precision of this method make it a valuable tool for the quantitative study of magneto-optic effects and for the characterization of optical materials in both educational and research settings.

Faraday rotation has various practical applications in science and technology. It is extensively used in telecommunications to rotate the polarization of light in optical fibers, aiding signal transmission and reducing interference. In astronomy, it helps measure magnetic fields in distant celestial objects. Faraday rotators are vital components in laser systems, enabling precise control of laser beam polarization. Additionally, it finds applications in magneto-optical data storage devices, where the rotation of light polarization is employed for data encoding. Faraday rotation plays a crucial role in studying magnetic materials, contributing to our understanding of their properties.

References

- [1] Michael Faraday. Experimental researches in electricity. nineteenth series. *Philosophical Transactions of the Royal Society of London*, 1846.
- [2] F. J. Loeffler. A faraday rotation experiment for undergraduate physics laboratory. *American Journal of Physics*, 51(7):623–625, 1983.
- [3] A. M. Portis. *Electromagnetic Fields*. John Wiley & Sons, Year.
- [4] D. A. Van Baak. Resonant faraday rotation as a probe of atomic dispersion. *American Journal of Physics*, 64(6):789–797, 1996.